## **A clinical overview of Roche SARS-CoV-2 antibody tests**

*Elecsys® Anti-SARS-CoV-2 (qualitative) Assay*

*Elecsys® Anti-SARS-CoV-2 S (semi-quantitative) Assay*


## Targeting high affinity antibodies is important for detecting previous exposure to SARS-CoV-2

When determining individual exposure and community prevalence of COVID-19, it is important to select tests that minimize the risk of false-positive results, which could infer a level of immune response or protection that is not actually present.

When COVID-19 prevalence is low, CDC guidance is to consider employing two independent tests in sequence when the first test yields a positive result in order to minimize false-positive test results and improve positive predictive value.<sup>1</sup> Additionally, the FDA suggests a second test using a different viral protein to make an informed decision on whether or not an individual has had a prior infection or truly has antibodies to the virus.<sup>2</sup>



The Roche **Elecsys® Anti-SARS-CoV-2 S (semi-quantitative)** and **Anti-SARS-CoV-2 (qualitative)** immunoassays detect high affinity antibodies against the spike and nucleocapsid proteins, respectively, to provide high specificity—delivering results you can trust to inform clinical decisions and support patient care.

### Elecsys® Anti-SARS-CoV-2 Assays: Performance by Design

Using a unique double-antigen sandwich assay format, Roche Elecsys® Anti-SARS-CoV-2 assays provide you with the flexibility to detect the relevant antibodies to achieve high specificity: high affinity, mature antibodies — predominantly immunoglobulin G (IgG), avoiding the drawbacks of immunoglobulin M (IgM)-only tests.


Factors contributing to high Negative Percent Agreement


Elecsys® Anti-SARS-CoV-2 S (semi-quantitative) Assay  
**99.98% Negative Percent Agreement**  
(N=5,991)<sup>6</sup>

Elecsys® Anti-SARS-CoV-2 (qualitative) Assay  
**99.80% Negative Percent Agreement**  
(N=10,453)<sup>7</sup>

### Estimated course of markers in SARS-CoV-2 infection<sup>8</sup>

It's clear that SARS-CoV-2 is no ordinary disease, and it doesn't behave like one either. While **IgM** antibodies often appear before **IgG** in other diseases, research shows that **IgG and IgM largely appear at the same time in SARS-CoV-2**. Since **IgM antibodies decline earlier than IgG antibodies**, are **less specific**, have a **high risk of false positives**, and are **not suitable to diagnose acute infection**, there is limited value in detecting class specific antibodies for SARS-CoV-2. This makes **IgG, total, or high affinity antibodies better suited for the use cases of SARS-CoV-2 antibody testing**.

## Test Principle

Elecsys® Anti-SARS-CoV-2 immunoassays utilize an in-solution double-antigen sandwich (DAGS) format, which requires the simultaneous binding of two differently labeled antigens via high affinity interactions to generate a signal. This reduces the probability of non-specific reactions with low affinity antibodies and cross-reaction with antibodies from other coronaviruses.

### Elecsys® Anti-SARS-CoV-2 S (semi-quantitative) Assay

Immunoassay for the **semi-quantitative** determination of antibodies against the SARS-CoV-2 **spike (S) protein**

#### Intended Use

Elecsys® Anti-SARS-CoV-2 S for use on the **cobas e** analyzers is an electrochemiluminescence immunoassay intended for the qualitative and semi-quantitative detection of antibodies to SARS CoV 2 spike (S) protein receptor binding domain (RBD) in human serum and plasma (lithium heparin, dipotassium-EDTA, tripotassium-EDTA, and sodium citrate). The Elecsys® Anti-SARS-CoV-2 S assay is intended as an aid in identifying individuals with an adaptive immune response to SARS-CoV-2, indicating recent or prior infection.<sup>6</sup>

### Elecsys® Anti-SARS-CoV-2 (qualitative) Assay

Immunoassay for the **qualitative** detection of antibodies against SARS-CoV-2 **nucleocapsid (N) protein**

#### Intended Use

Elecsys® Anti-SARS-CoV-2 for use on the **cobas e** analyzers is an electrochemiluminescence immunoassay intended for the qualitative detection of antibodies to SARS-CoV-2 in human serum and plasma. The assay uses a recombinant protein representing the nucleocapsid (N) antigen for the determination of antibodies against SARS-CoV-2. The test is intended for use as an aid in identifying individuals with an adaptive immune response to SARS-CoV-2, indicating recent or prior infection.<sup>7</sup>

## Why use Elecsys® SARS-CoV-2 antibody tests?

### Measure adaptive immune response

- Indicate recent or prior infection

### Determine seroprevalence in a given population in high, moderate and low prevalence settings

### Determine seroprevalence in a given population in very low to moderate seroprevalence settings

- Perform sequential testing to improve specificity in low prevalence settings

### Increase sensitivity by combining Anti-S and Anti-N

- Identify S-only AND N-only responders

### Determine antibody levels in convalescent plasma donations


- Assist with preparation of therapeutic Ig formulations

### Evaluate vaccine-induced immune response\*

- Support vaccine evaluation studies
- Semi-quantitatively measure the immune response to RBD-containing vaccines
- Determine pre-vaccination immune status in order to administer vaccine to sero-naïve individuals

#### Elecsys® Anti-SARS-CoV-2 S (semi-quantitative) Assay

#### Elecsys® Anti-SARS-CoV-2 (qualitative) Assay

\* Subject to vaccine emergency use authorization and CDC guidelines.

Assay Characteristics	Elecsys® Anti-SARS-CoV-2 S Assay	Elecsys® Anti-SARS-CoV-2 Assay
<b>Specimen type</b>	Serum collected using standard sampling tubes or tubes containing separating gel; Li-heparin, K2-EDTA-, K3-EDTA-, and sodium citrate plasma; plasma tubes containing separating gel can be used.	Serum collected using standard sampling tubes or tubes containing separating gel. Li-heparin, K-EDTA and K-EDTA plasma as well as plasma tubes containing separating gel.
<b>Sample volume</b>	12 to 20 µL, depending on system used	
<b>Testing time</b>	18 minutes	
<b>Traceability</b>	Internal Roche standard for anti-SARS-CoV-2-S consisting of monoclonal antibodies. 1 nM of these antibodies correspond to 20 U/mL of the Elecsys® Anti-SARS-CoV-2 S assay.	Internal positive and negative samples
<b>Linear range</b>	0.4 to 250 U/mL	(qualitative)
<b>Result interpretation</b>	<0.8 U/mL = negative, ≥0.8 U/mL = positive	COI (cutoff index) <1.0 = non-reactive, COI ≥1.0 = reactive
<b>Positive Percent Agreement (PPA)</b>	<b>1,485 samples</b> from 331 symptomatic patients with a PCR confirmed SARS-CoV-2 infection were tested. 233 samples ≥15 days after diagnosis with PCR. 225 of these 233 samples ≥0.8 U/mL considered positive, PPA of <b>96.6%</b> (95% confidence interval: 93.35 – 98.51%) in this sample cohort.	<b>496 samples</b> from 102 symptomatic patients with a PCR confirmed SARS-CoV-2 infection were tested. One or more sequential specimens were collected after PCR confirmation at various time points. PPA of <b>99.5%</b> (95% confidence interval: 97.0 – 100%) in this sample cohort.
<b>Negative Percent Agreement (NPA)</b>	<b>5,991 samples</b> drawn before Oct. 2019 were tested. Overall NPA in this cohort of potentially cross-reactive samples was <b>99.98%</b> (95% confidence interval: 99.91 – 100%).	<b>10,453 samples</b> drawn before Dec. 2019 were tested. Overall NPA in this cohort of potentially cross-reactive samples was <b>99.80%</b> (95% confidence interval: 99.69 – 99.88%).
<b>Cross-reactivity (analytical specificity)</b>	Showed <b>no cross-reactivity</b> in <b>1,100</b> potentially cross-reactive samples including 7 samples from patients with MERS-CoV, 21 samples from patients with symptoms of common cold, 94 samples from other coronaviruses, and 978 pre-pandemic samples with reactivity for various other indications.	Out of <b>792</b> potentially cross-reacting samples, <b>4 samples</b> showed reactivity in the Elecsys® Anti-SARS-CoV-2 assay resulting in an overall specificity in this cohort of <b>99.5%</b> (95% confidence interval: 98.63–99.85 %).

## REFERENCES

1. <https://www.cdc.gov/coronavirus/2019-ncov/lab/resources/antibody-tests-guidelines.html>. Accessed 12/08/2020.
2. <https://www.fda.gov/medical-devices/coronavirus-disease-2019-covid-19-emergency-use-authorizations-medical-devices/eua-authorized-serology-test-performance>. Accessed 12/08/2020.
3. Roche Diagnostics. SARS-CoV-2 antibody specificity white paper. June 29, 2020.
4. Callaway E. The race for coronavirus vaccines: a graphical guide. *Nature*. 2020;580:576–7.
5. Zhou G, Zhao Q. Perspectives on therapeutic neutralizing antibodies against the Novel Coronavirus SARS-CoV-2. *Int J Biol Sci*. 2020;16:1718–23.
6. Elecsys® Anti-SARS-CoV-2 S Assay. Package Insert 2020-12, V1.0; Material Numbers 09289267190 and 09289275190. Roche Diagnostics.
7. Elecsys® Anti-SARS-CoV-2 Assay. Package Insert 2020-07, V3.0; Material Numbers 09203095190 and 092030. Roche Diagnostics.
8. Adapted from Sethuraman et al. Estimated time intervals and rates of viral detection are based on data from several published reports. Because of variability in values among studies, estimated time intervals should be considered approximations and the probability of detection of SARS-CoV-2 infection is presented qualitatively. Sethuraman N, Jeremiah SS, Ryo A. Interpreting diagnostic tests for SARS-CoV-2. *JAMA*. Published online May 06, 2020. doi:10.1001/jama.2020.8259.
9. CDC, Interim Guidelines for COVID-19 Antibody Testing: <https://www.cdc.gov/coronavirus/2019-ncov/lab/resources/antibody-tests-guidelines.html>.
10. Deeks JJ, et al. Antibody tests for identification of current and past infection with SARS-CoV-2 (Review). *Cochrane Database of Systematic Reviews* 2020, Issue 6. Art. No.: CD013652. doi: 10.1002/14651858.CD013652.
11. Cheng MP, et al. Serodiagnostics for severe acute respiratory syndrome-related coronavirus-2: a narrative review. *Ann Intern Med*. 2020. doi:10.7326/M20-2854.
12. Sethuraman N, et al. Interpreting diagnostic tests for SARS-CoV-2. *JAMA*. 2020;323:2249–51.
13. To K, et al. Temporal profiles of viral load in posterior oropharyngeal saliva samples and serum antibody responses during infections by SARS-CoV-2: an observational cohort study. *The Lancet*. 2020;20(5):P565–74. [https://doi.org/10.1016/S1473-3099\(20\)30196-1](https://doi.org/10.1016/S1473-3099(20)30196-1).
14. Long Q, et al. Antibody responses to SARS-CoV-2 in COVID-19 patients: the perspective application of serological tests in clinical practice. medRxiv 2020. <https://www.medrxiv.org/content/10.1101/2020.03.18.20038018v1>.
15. Xiang F, et al. Antibody detection and dynamic characteristics in patients with coronavirus disease 2019. *Clin Infect Dis*. 2020;71(8):1930–34. <https://doi.org/10.1093/cid/ciaa461>.
16. Hou H, et al., Detection of IgM and IgG antibodies in patients with coronavirus disease 2019. *Clin Transl Immunology*. 2020;9(5): e01136. <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC7202656/pdf/>

Not for screening of donated blood. This test has not been FDA cleared or approved. This test has been authorized by FDA under an EUA for use by authorized laboratories. This test has been authorized only for detecting the presence of antibodies against SARS-CoV-2, not for any other viruses or pathogens. This test is only authorized for the duration of the declaration that circumstances exist justifying the authorization of the emergency use of in vitro diagnostic tests for detection and/or diagnosis of COVID-19 under section 564(b)(1) of the Act, 21 U.S.C. § 360bbb- 3(b)(1), unless the authorization is terminated or revoked sooner.

**For more information, contact your sales representative  
or visit [go.roche.com/cobas-SARS-CoV-2](https://go.roche.com/cobas-SARS-CoV-2)**

COBAS, COBAS E, and ELECSYS are trademarks of Roche.

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#### Referencias

1. <https://www.cdc.gov/coronavirus/2019-ncov/lab/resources/antibody-tests-guidelines.html>. Accessed 12/08/2020.
  2. <https://www.fda.gov/medical-devices/coronavirus-disease-2019-covid-19-emergency-use-authorizations-medical-devices/eua-authorized-serology-test-performance>. Accessed 12/08/2020.
  3. Roche Diagnostics. SARS-CoV-2 antibody specificity white paper. June 29, 2020.
  4. Callaway E. The race for coronavirus vaccines: a graphical guide. *Nature*. 2020;580:576-7.
  5. Zhou G, Zhao Q. Perspectives on therapeutic neutralizing antibodies against the Novel Coronavirus SARS-CoV-2. *Int J Biol Sci*. 2020;16:1718-23.
  6. Elecsys® Anti-SARS-CoV-2 S Assay. Package Insert 2020-12, V1.0; Material Numbers 09289267190 and 09289275190. Roche Diagnostics.
  7. Elecsys® Anti-SARS-CoV-2 Assay. Package Insert 2020-07, V3.0; Material Numbers 09203095190 and 092030. Roche Diagnostics.
  8. Adapted from Sethuraman et al. Estimated time intervals and rates of viral detection are based on data from several published reports. Because of variability in values among studies, estimated time intervals should be considered approximations and the probability of detection of SARS-CoV-2 infection is presented qualitatively. Sethuraman N, Jeremiah SS, Ryo A. Interpreting diagnostic tests for SARS-CoV-2. *JAMA*. Published online May 06, 2020. doi:10.1001/jama.2020.8259.
  9. CDC, Interim Guidelines for COVID-19 Antibody Testing: <https://www.cdc.gov/coronavirus/2019-ncov/lab/resources/antibody-tests-guidelines.html>.
  10. Deeks JJ, et al. Antibody tests for identification of current and past infection with SARS-CoV-2 (Review). *Cochrane Database of Systematic Reviews* 2020, Issue 6. Art. No.: CD013652. doi: 10.1002/14651858.CD013652.
  11. Cheng MP, et al. Serodiagnostics for severe acute respiratory syndrome-related coronavirus-2: a narrative review. *Ann Intern Med*. 2020. doi:10.7326/M20-2854.
  12. Sethuraman N, et al. Interpreting diagnostic tests for SARS-CoV-2. *JAMA*. 2020;323:2249–51.
  13. To K, et al. Temporal profiles of viral load in posterior oropharyngeal saliva samples and serum antibody responses during infections by SARS-CoV-2: an observational cohort study. *The Lancet*. 2020;20(5):P565-74. [https://doi.org/10.1016/S1473-3099\(20\)30196-1](https://doi.org/10.1016/S1473-3099(20)30196-1).
  14. Long Q, et al. Antibody responses to SARS-CoV-2 in COVID-19 patients: the perspective application of serological tests in clinical practice. medRxiv 2020. <https://www.medrxiv.org/content/10.1101/2020.03.18.20038018v1>.
  15. Xiang F, et al. Antibody detection and dynamic characteristics in patients with coronavirus disease 2019. *Clin Infect Dis*. 2020;71(8):1930-34. <https://doi.org/10.1093/cid/ciaa461>.
- Hou H, et al., Detection of IgM and IgG antibodies in patients with coronavirus disease 2019. *Clin Transl Immunology*. 2020;9(5): e01136. <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC7202656/pdf/>

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© 2020 Roche. MC-US-08237-1220
Roche Diagnostics Corporation